BioGenex8810 Kato Road, Suite 100E & 200E,
Fremont, CA94538Tel : +1(800)421-4149,
Fax: +1(510)824-1490,
support@biogenex.com

ECIREP

EmergoEurope
Prinsessegracht
20 2514 AP
The Hague,,
The Netherlands

COUNTER STAINS

Doc. No. 932-HK100E-4 Rev. No. 1
Release Date: 10-Aug-2020

ENGLISH

Reagents Provided: One of the following

Catalog No.	Counter Stain
HK100-5K, HK100-9K	Hematoxylin

Intended Use

BioGenex Counter stains are intended for in vitro diagnostic use to provide color contrast in formalin-fixed tissue sections in an immunohistochemistry (IHC) or in situ hybridization (ISH) assay.

Summary and Explanation

Counter stains provide complementary information regarding the overall morphology or histology of the cell or tissue to that provided by the primary stain. Thoughtful selection of a combination of the primary stain and counter stain can result in significant enhancement of staining signal. BioGenex provides several effective counter stain options that have a contrasting and subtle color, which do not interfere with the primary stain.

Principles of the Procedure

The counter stains take advantage of the chemical nature of the cellular components and react with them. The stains may be specific for either the acidic or the basic components of cells.

Reagents Required but Not Supplied

All reagents required for counter staining are provided. See antibody and detection kit datasheets for complete set of reagents required for immunohistochemistry (IHC) or in situ hybridization (ISH) procedures.

Storage and Handling

Store all reagents at 2-8°C. Do not use after expiration dates as indicated on the reagent labels.

Specimen Preparation

Paraffin embedded and formalin fixed tissue sections 4-6 micron thick are suitable for use with enzyme pretreatment.

Precautions

HK100-XX:

Xi;R36/38 = Irritating to eyes and skin.

S26 S37 S24/25 S60 = In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable gloves. Avoid contact with skin and eyes. This material and its container must be disposed of as hazardous waste.

Staining Protocol

The following table provides the recommended treatments with counter stain. Refer to the antibody and detection kit datasheets for complete staining protocols. User should validate any other conditions. For automated use follow instructions in operator's manual.

Counter stain	BioGenex Recommendations

Hematoxylin	Stain slides for 1 to 5 minutes, rinse with distilled water. Dehydrate sections in a series of graded alcohols and clear in xylene or xylene substitute, then mount. Compatible with AEC, DAB and Fast Red chromogens.
-------------	--

Quality Control

Refer to the appropriate detection system package inserts for guidance on general quality control procedures.

Trouble shooting

- Run appropriate positive and negative controls to ensure the optimal signal.
- Refer to the troubleshooting section in the package inserts of BioGenex Super Sensitive Detection Systems (or other equivalent detection systems) for remedial actions on detection system related issues, or contact BioGenex Technical Service Department at 1 800-421-4149 or customer.service@biogenex.com to report unusual staining.

Expected Results

Counter stain	Expected Results
Hematoxylin	Nucleus stains bluish purple.
Methyl Green	Nucleus stains light green
Nuclear Fast Red	Nucleus stains red, cytoplasm stains pink.

Limitations of the procedure

Methyl green leaches in aqueous mounting medium so use only permanent mounting medium with methyl green counterstain. Improper dewaxing will lead to negative or weak signal.

Performance Characteristics

BioGenex has conducted studies to evaluate the performance of all its counter stains using several BioGenex antibodies and detection systems. The counter stains have been found to show color contrasts with respect to the primary stain thereby enhancing the staining pattern. BioGenex counter stains have shown reproducible and consistent results when used within a single run, between runs, between lots and wherever applicable between manual and automated runs. The products have been determined to be stable for the periods specified on the labels either by standard real time or accelerated testing methods. BioGenex ensures product quality through 100% quality control for all products released and through surveillance programs.

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Bibliography

- Sheehan DC and Hrapchak BB. Theory and Practice of Histotechnology. St. Louis: C.V. Mosby Co. 1980.
- Kiernan JA. Histological and Histochemical Methods: Theory and Practice. New York: Pergamon Press 1981.
- Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. Lab Med 1983; 14:767-770.
- Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immuno histochemistry. Am J ClinPathol 1980 May;73(5):626-632.

	Representative in the European Community MandatarionellaComunitàEuropea Bevollmächtigter in der EuropäischenGemeinschaft Representantautorizadoen la ComunidadEuropea
	In Vitro Diagnostic Medical Device Dispositivo medico-diagnostico in vitro In Vitro Diagnostikum Productsanitario para diagnóstico in vitro
	Cons ult Instructions for use Consultarele istruzioni per l'uso Gebrauchsanweisungbeachten Consulte las instruccionesde uso
	Temperature Limitation Limiti di temperatura ZulässigerTemperaturbereich Limite de temperatura
	Manufacturer Fabricante Hersteller Fabricante