

DATA SHEET

Anti-CMV (Cytomegalovirus) Mouse Monoclonal Antibody

Catalog No.	Description
AM254-5M	6 ml of Prediluted Antibody
MU254-UC	1 ml of Concentrated Antibody
MU254-5UC	0.5 ml of Concentrated Antibody

Analyte Specific Reagent. Analytical and performance characteristics are not established.

Clone

BM204

Immunoglobulin Class

Mouse IgG1

Specifications

This antibody stains Cytomegalovirus in formalin-fixed, paraffin-embedded tissue sections by immunohistochemical techniques.

Storage

Store at 2-8°; do NOT freeze. Do not use after expiration date on vial.

Source and Format

Mouse Monoclonal Antibody to Cytomegalovirus from cell culture supernatant, diluted in PBS, pH 7.6, containing 1% BSA and 0.09% sodium azide.

Precautions

For Professional use. This product contains no hazardous material at a reportable concentration according to U.S. 29 CFR 1910.1200, OSHA Hazard Communication Standard and EC Directive 91/155/EC. However, this product contains sodium azide, at concentrations of less than 0.1%. Sodium azide is not classified as a hazardous chemical at the product concentrations. However, toxicity information regarding sodium azide at product concentrations has not been thoroughly investigated. Sodium azide may react with lead or copper plumbing to form highly explosive metal azides. Upon

disposal, flush with large volumes of water to prevent azide build-up in plumbing (Center for Disease Control, 1976, National Institute for Occupational Safety and Health, 1976). For more information, a Safety Data Sheet (SDS) for sodium azide in pure form is available upon request. Do not pipette reagents by mouth, and avoid contact of reagents and specimens with skin and mucous membranes. If reagents or specimens come in contact with sensitive area, wash with copious amounts of water. Minimize microbial contamination of reagents or else increase in nonspecific staining may occur.

Quality Control

Each lot of this antibody is tested by immunohistochemistry for Quality Control purposes. Refer to the BioGenex Quality Control Testing Conditions table for additional information.

References

Ruebner B H, et al. Am J Pathol 48:971-989, 1966.
Manual of Clinical Microbiology, Balows Am et al., eds. American Society of Microbiology, Washington, D.C., 1991:829-837.
Stern H and Tucker S.M., Br Med J 2:268, 1973.
Monif GR, et al. J. Pediatr 80:17, 1972.
Center for Disease Control. Center for Disease Control Manual Guide- Safety Managements, No. CDC-22, Atlanta, Georgia. April 30,1976.
Natal Institute for Occupational Safety and Health, (NIOSH), Rockville, MD. Publication No. 78-127, 1976.

BioGenex Quality Control Testing Conditions

Parameter	Conditions Used
Control Tissue	CMV INF. LUNG
Tissue Type	Formalin-Fixed, Paraffin-Embedded tissue sections (FFPE)