BioGenex
48810 Kato Road, Suite 100E & 200E,
Fremont, CA 94538

CE
EC/REP
Emergo Europe,
Prinsessegracht
20 2514 AP
The Hague,
The Netherlands

Tel : +1 (800) 421-4149,
Fax: +1 (510) 824-1490,
support@biogenex.com

LINKS AND LABELS

Doc. No. 932-HK268E-4 Rev. No. I
Release Date: 10-Aug-2020

ENGLISH

Intended Use

Biotinylated secondary antibodies (Links) and streptavidin-enzyme conjugates (Labels) are intended for in vitro diagnostic use for the detection of primary antibodies in paraffin-embedded tissue sections, cryostat sections, cell monolayers and whole mounts.


Summary and Explanation

Ready to use links and labels are optimized for use with ready-to-use (RTU) antibodies. The links are available as either species specific or Multilink® secondary antibodies. Multilink® permits detection of mouse IgG and IgM and rabbit immunoglobulins, providing a universal link reagent for use with most monoclonal and polyclonal antibodies. Concentrated links and labels are applied with recommended dilutions for use with antibodies. The working solutions are stable for 6 months when diluted and stored under recommended conditions. The concentrated links should be diluted with Link Diluent (Cat. No. HK165-5K, except for concentrated goat link which should be diluted with HK156-5K). Concentrated peroxidase labels should be diluted in streptavidin peroxidase Diluent (Cat. No. HK157-5K), while concentrated alkaline phosphatase labels should be diluted with Common Antibody Diluent (Cat. No. HK156-5K)

Principles of the Procedure

Antigens in tissues and cells are detected by a two-stage process—the binding of the primary antibody to a specific epitope, and the subsequent detection of this binding by a colorimetric reaction. Tissues or cell preparations are frozen or fixed, sectioned, and attached to slides. The sections are then dewaxed if paraffin-embedded, treated with an antigen retrieval solution if required, may be blocked with required blocking solutions and then incubated with a primary antibody (1). The bound primary antibody is detected by the serial addition of biotinylated secondary antibody or Link (2), Streptavidin-enzyme conjugate or Label (3), and substrate (4). When adequate color development is seen, the slides are washed in water to stop the reaction, counterstained, and covered with a mounting medium. With this detection system, signal amplification is achieved by the binding of multiple units of secondary antibody to each primary antibody molecule, followed by the binding of multiple enzyme conjugates to each secondary antibody, and finally the enzymatic conversion of the substrate.

Reagents Provided

One of the following:

1 vial (6, 17, 50, or 100 ml) of ready-to-use Link (biotinylated anti-immunoglobulin in PBS with carrier protein, and 0.09% sodium azide).

1 vial (6, 17, 50, or 100 ml) of ready-to-use HRP Label (horseradish peroxidase-conjugated streptavidin in PBS with carrier protein and 0.1% Proclin 300).

1 vial (6, 17, 50 or 100 ml) of ready-to-use AP Label (alkaline phosphatase-conjugated streptavidin in PBS with carrier protein and 0.09% sodium azide).

1 vial (5 ml) of Concentrated Link (biotinylated anti-immunoglobulin in PBS with carrier protein, 0.09% sodium azide).

1 vial (5 ml) of Concentrated HRP Label (horseradish peroxidase-conjugated streptavidin in PBS with carrier protein and 0.1% Proclin 300).

1 vial (5 ml) of Concentrated AP Label (alkaline phosphatase-conjugated streptavidin in PBS with carrier protein and 0.09% sodium azide).

Cat. No.	Description
HK335-5ME	Super Sensitive Mouse Link (biotinylated goat anti-mouse Ig)
HK336-5R, 9R	Super Sensitive Rabbit Link (biotinylated goat anti-rabbit Ig)
HK340-5K, 9K	Super Sensitive Multilink® (biotinylated goat anti-Ig; for use with mouse and rabbit primary antibodies)
HK330-5K, 9K	Super Sensitive HRP Label (peroxidase-conjugated streptavidin)
HK331-5K, 9K	Super Sensitive AP Label (alkaline phosphatase-conjugated streptavidin)
HK268-UKE	Concentrated Multilink® (biotinylated anti-Ig; for use with mouse, rabbit, guinea pig, and rat primary antibodies)
HK320-UKE	Concentrated HRP Label (peroxidase-conjugated streptavidin)
HK321-UKE	Concentrated AP Label (alkaline phosphatase-conjugated streptavidin)

For Use in Manual Procedures and with BioGenex Automated Staining Systems

Reagents Required but Not Supplied

EZ-DeWax™ Solution (optional), Antigen retrieval solutions (optional), Blocking solutions, Primary antibodies, Link/ Label, substrates, Counterstain and mounting medium.

Storage and Handling

Store all reagents at 2-8°C. Do not use after expiration dates as indicated on the reagent labels.

Specimen Preparation

Tissues fixed in 10% (v/v) formalin prior to paraffin embedding are suitable for use. Consult references (Kiernan, 1981; Sheehan & Hrapchak, 1980) for further details on specimen preparation.

Precautions

Links and alkaline phosphatase Labels contain sodium azide at concentrations of less than 0.1%. Sodium azide is not classified as a hazardous chemical at the concentration of this product. However, toxicity information regarding sodium azide at the product’s concentration has not been thoroughly investigated. For more information, a Material Safety Data Sheet (MSDS) for sodium azide in pure form is available upon request.

Buffers containing sodium azide should not be used to dilute peroxidase labels. Links and alkaline phosphatase labels may be diluted in Common Antibody Diluent, and peroxidase labels may be diluted in Streptavidin Peroxidase Diluent, both available from BioGenex. These reagents can also be titrated to different sensitivity levels to match the user’s application; however, this and use of other diluents should be tested extensively and validated in the user’s laboratory.

Preparation of Working Solutions

Dilute Concentrated Links and Labels 1:20 with the appropriate recommended diluent.

Staining Protocol

- BioGenex ready-to-use Links and Labels have been optimized for use with BioGenex ready-to-use primary antibodies.
- After application of the primary antibody, tissue sections should be rinsed well with PBS at room temperature prior to addition of the Link.
- Carefully wipe excess liquid from around the tissue section. Apply sufficient Link (1-2 drops) to cover the section entirely, and incubate for 20 minutes at room temperature or 5 minutes at 37°C. Rinse well with PBS.
- Carefully wipe excess liquid from around the tissue section. Apply sufficient Label (1-2 drops) to cover the section entirely, and incubate for 20 minutes at room temperature or 5 minutes at 37°C. Rinse well with PBS.
- Sections are now ready for application of a suitable enzyme substrate.
- Note: refer to your BioGenex Automated Stainer operator’s manual for use of these products in automated procedures.

Quality Control

Refer to the appropriate detection system package inserts for guidance on general quality control procedures.

Trouble shooting

Run appropriate positive and negative controls to ensure the optimal signal.

Refer to the troubleshooting section in the package inserts of BioGenex Super Sensitive Detection Systems (or other equivalent detection systems) for remedial actions on detection system related issues, or contact BioGenex Technical Service Department at (925) 275-0550 to report unusual staining

Expected Results

The colored pigment deposition will be seen at the site of antigen-antibody reaction. The color signal can be visualized using light microscopy.

Limitations of the procedure

Immunohistochemistry (IHC) is a complex technique involving both histological and immunological detection methods. Tissue processing and handling prior to immunostaining can also cause inconsistent results. Variations in fixation and embedding or the inherent nature of the tissue may cause variations in results (Nadji and Morales, 1983). Endogenous peroxidase activity or pseudoperoxidase activity in erythrocytes and endogenous biotin may cause non-specific staining depending on detection system used. Tissues containing hepatitis B surface antigen (HBsAg) may give false positive with horse radish peroxidase systems. (Omata M et al, 1980). Improper counterstaining and mounting may compromise the interpretation of results.

Performance Characteristics

BioGenex has conducted studies to evaluate the performance of all its links and labels using several BioGenex antibodies and detection systems. BioGenex links and labels have shown reproducible and consistent results when used within a single run, between runs, between lots and wherever applicable between manual and automated runs. The products have been determined to be stable for the periods specified on the labels either by standard real time or accelerated testing methods. BioGenex ensures product quality through 100% quality control for all products released and through surveillance programs.

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Bibliography

- Sheehan DC and Hrapchak BB. Theory and Practice of Histotechnology. St. Louis: C.V. Mosby Co. 1980.
- Kiernan JA. Histological and Histochemical Methods: Theory and Practice. New York: Pergamon Press 1981.
- Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. Lab Med 1983; 14:767-770.
- Omata M, Liew CT, Ashcavi M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen. A possible source of error in immuno histochemistry. Am J Clin Pathol 1980 May;73(5):626-632.

EC REP	Representative in the European Community <p>Mandatario nella Comunità Europea</p> Bevollmächtigter in der Europäischen Gemeinschaft <p>Representante autorizado en la Comunità Europea</p>
IVD	In Vitro Diagnostic Medical Device <p>Dispositivo medico-diagnostico in vitro</p> In Vitro Diagnostikum <p>Producto sanitario para diagnóstico in vitro</p>
i	Consult Instructions for use <p>Consultare le istruzioni per l'uso</p> Gebrauchsanweisung beachten <p>Consulte las instrucciones de uso</p>
T	Temperature Limitation <p>Limiti di temperatura</p> Zulässiger Temperaturbereich <p>Limite de temperatura</p>
M	Manufacturer <p>Fabbricante</p> Hersteller <p>Fabri cante</p>