

## DATA SHEET

### Hsa-miR-375Probe

<b>Catalog No.</b> <b>HM375-100</b>	<b>Description</b> One vial of 0.650 ml of probe in hybridization buffer
--	---

**Analyte Specific Reagent. Analytical and performance characteristics are not established.**

---

Doc. No. 932-HM375-100      Rev : D  
Date of release: 17-Aug-2020

#### Description

The Hsa-miR-375 probe has been designed from mature human miR-375 sequence. This fluoresceinated probe is provided in a hybridization buffer for localization of miRNA in FFPE tissue by *In Situ* hybridization.

#### Specifications

The Hsa-miR-375 identifies mature miR-375 sequences in formalin-fixed, paraffin-embedded human tissues and/or freshly prepared frozen tissues by *in situ* hybridization. This probe does not react with normal human mRNA or nuclear DNA present in tissues.

#### Storage and Handling

Store the reagent at 2-8 °C. Do not freeze. Do not use the reagent after expiration date on vial. The reagent must be brought to room temperature before use. (Important! The presence of precipitates induces background staining).

#### Precautions:

For professional use the probe contains formamide. Formamide is classified as a teratogen. Pregnant workers should keep exposure to a minimum. Avoid inhalation, ingestion, and contact with unprotected skin. If skin contact occurs, wash thoroughly with soap and water. For more information, refer to the Material Safety Data Sheet, which is available upon request.

#### Quality Control

Each lot of this micro RNA probe is tested by *In Situ* hybridization for Quality Control purposes. Refer to the BioGenex Quality Control Testing Conditions table for additional information.

#### References

1. Lorio MV and Croce CM. (2012). MicroRNA dysregulation in cancer: diagnostics, monitoring and therapeutics. A comprehensive review. **EMBO Mol Med** 4, 143–159.
2. Chen PS, Su JL, and Hung MC. (2012). Dysregulation of Micro RNAs in cancer. **Journal of Biomedical Science**, 19:90.
3. Nuovo GJ. (2008). In situ detection of precursor and mature microRNAs in paraffin embedded, formalin fixed tissues and cell preparations. **Methods** 44,39–46.

# **BioGenex**

48810 Kato Road, Suite 100E & 200E  
Fremont, CA 94538  
Tel : +1 (800) 421-4149, Fax: +1 (510) 824-1490,  
support@biogenex.com

4. Song R. et al. (2010). In situ hybridization detection of microRNAs. **Methods Mol Biol.** 629, 287-94.
5. Hartz, P. A. Personal Communication. Baltimore, Md. 8/20/2010.
6. Poy, MN., Hausser, J, Trajkovski M, Braun M, Collins S., Rorsman, P., Zavolan, M., Stoffel, M. (2009). miR-375 maintains normal pancreatic alpha- and beta-cell mass. **Proc. Nat. Acad. Sci.** 106, 5813-5818.
7. Zhang, X., et al. (2011). Combination of hsa-miR-375 and hsa-miR-142-5p as a predictor for recurrence risk in gastric cancer patients following surgical resection. **Annals of oncology** 22 (10), 2257-2266.

## **BioGenex Quality Control Testing Conditions**

<b>Parameter</b>	<b>Conditions used</b>
Control Tissue	COLON, HEMANGIOMA. KIDNEY (FB-HM375)
Tissue Type	Formalin-fixed, paraffin-embedded cancer tissues