



49026 Milmont Drive, Fremont, CA 94538
Tel : +1 (800) 421-4149, Fax: +1 (510) 824-1490,
support@biogenex.com

DATA SHEET

***i6000*[™] and OptiMax[®] Pipette Tips**

Doc. No.932-XT011, Rev. No. C
Date of Release: 06-October, 2011

Intended Use

This product is intended for laboratory use only. BioGenex pipette tips are optimized for use on BioGenex OptiMax[®] or *i6000*[™] Staining Systems.

Specifications

BioGenex pipette tips are made of high-quality polypropylene and are RNase- and heavy metals-free when untampered. Inner surface is extremely smooth and requires minimum wetting. Each tip is capable of dispensing 1ml (1000) maximum volume.

Items Supplied

One of the following:

Cat. No.	Description
XT011-PT	1 case (1000 tips) of 1ml Pipette tips for OptiMax [®]
XT100-PT	1 box (100 tips) of 1ml Pipette tips for OptiMax [®]
XT104-05X	1 case (960 tips) of 1ml Pipette tips <i>i6000</i> [™]
XT105-01X	1 box (192 tips) of 1ml Pipette tips for <i>i6000</i> [™]
XT106-PT	1 case (1000 tips) of 1ml Pipette tips for OptiMax [®] Upgrade
XT107-PT	1 box (100 tips) of 1ml Pipette tips for OptiMax [®] Upgrade

Note:

1. Do not touch pipette tips with bare hands to avoid RNase contamination.
2. Pipette tip box should be properly aligned on the instrument. Improper alignment may lead to display of "Error" messages.